

Basic Sizing Analysis

Quick Reference Card


This quick reference card describes an analysis workflow for Peak Scanner™ Software Version 1.0 using data included with the software installation.

Note: This workflow does not include instructions for editing or managing size standards or analysis methods.

Launch the Peak Scanner Software

- Click  (Peak Scanner™ Software)
- or
- Select Start ▶ Programs ▶ Applied Biosystems ▶ Peak Scanner v1.0 ▶ Peak Scanner v1.0.

Components to the Peak Scanner Software Main Window



A	Application Toolbar
B	Task Pane
C	Page
D	Project View

E	Task Action Toolbar
F	Panel Action Toolbar
G	Panel

Peak Scanner™ Software Workflow

Launch the Peak Scanner Software

Select a new or existing project

Add sample files

Select a size standard and an analysis method


Click Analyze

Check plot view and sizing view panels to evaluate the size quality

Export or print the data

Select a New or Existing Project

Create a new project in the left pane or open a project in the right pane of the Start View window.



To create a new project:


1. Click Start New Project.
2. Add samples.
3. Assign analysis methods and size standards to each sample.
4. Analyze the samples to view data and sizing information.

To open a project:

Select from the list of displayed projects in your Application data folder or browse to the project location.

Add Sample Files

1. Select Add Files in the panel action toolbar.
2. Browse to the .fsa files you want to add to the project.
3. Select either the folder containing the files or the individual files to be added.
4. Click Add selected files, then click OK in the Add Sample Files dialog box.



You can also add files to a project by dragging the files/folders directly into the open application.

Select Size Standard and Analysis Method

Select a Size Standard

Under the column heading Size Standard, select the drop-down list to choose the size standard used for each sample.

Add Files Remove Files Size Standard						
Size Standard	Analysis Method	SQI	Offscale	Quality	UD1	UD
GS600LIZ	Sizing Default NPP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
GS600LIZ_3730	Sizing Default	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		

Select an Analysis Method

Under the column heading Analysis Parameter, select the drop-down menu to choose the analysis method.


Add Files Remove Files Size Standard						
Size Standard	Analysis Method	SQI	Offscale	Quality	UD1	UD
GS600LIZ	Sizing Default	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
GS600LIZ	Sizing Default NPP Sizing Default My AM	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		

The two default analysis methods are provided with the software:

- Sizing Default —Use when primers have not been removed from the sample.
- Sizing Default NPP (NPP= no primer peak(s))— Used when samples do not contain primers.


Click Analyze

- In the Task Pane, click Analyze.



or

- In the Application toolbar, click Analysis
► Analyze All.




Check Plot View and Sizing View Panels to Evaluate Size Quality

Check Plot View Panel

Review electropherograms and annotate data using the Plot View Panel. The Plot View Panel is activated immediately after analysis of sample files but you can access it using the task action toolbar from multiple pages.

After data analysis, launch the plot panel in a new window to maximize the view. Minimize the Plot View Panel to gain more screen space for the Samples and Sizing Table View Panels.



Sizing Table

View peak information for the selected samples within the Samples View panel using the Sizing Table View panel. The Sizing table contains a list of all peaks for the selected samples from the Samples Table with properties that are associated with each peak.

Dye Color	Dye/Sample Peak	Sample File Name	Size	Height	Area in BP	Data Point	User Comments
49	Y_11	C006_E109_Run_C006_2002-05-10...	44,3763	522	400	2068	
50	Y_12	C006_E109_Run_C006_2002-05-10...	118,9394	488	300	3100	
51	Y_13	C006_E109_Run_C006_2002-05-10...	120,9395	352	166	3125	
52	Y_14	C006_E109_Run_C006_2002-05-10...	170,7026	292	179	3800	
53	Y_15	C006_E109_Run_C006_2002-05-10...	215,7185	601	433	4246	
54	Y_16	C006_E109_Run_C006_2002-05-10...	222,3916	460	343	4328	
55	Y_17	C006_E109_Run_C006_2002-05-10...	312,2798	200	167	5327	
56	Y_18	C006_E109_Run_C006_2002-05-10...	314,1392	319	269	5347	
57	R_1	C006_E109_Run_C006_2002-05-10...	12,5661	2804	2230	1587	
58	R_2	C006_E109_Run_C006_2002-05-10...	14,0212	2873	3617	1609	
59	R_3	C006_E109_Run_C006_2002-05-10...	14,6825	2317	1604	1619	
60	R_4	C006_E109_Run_C006_2002-05-10...	16,2037	3667	1473	1642	
61	R_5	C006_E109_Run_C006_2002-05-10...	17,1958	310	295	1657	
62	R_6	C006_E109_Run_C006_2002-05-10...	18,6508	318	298	1679	
63	R_7	C006_E109_Run_C006_2002-05-10...	19,5767	450	79	1693	

Export or Print the Data

Export

Export the currently selected data table or a combined table as either a:

- Tab Delimited format (*.txt)
or
- Comma separated values file (*.csv).

To export the Sizing or the Sample table:

- Select the panel containing the table to be exported, then click Export () from the application toolbar.
- Select Export.

Note: The Export selected table option is enabled only when viewing the Samples or the Sizing Table View panels.

The Combined table contains information displayed in both the Sample and Sizing tables.

To export the Combined table:

- Click Export () in the application toolbar.
- Click Export Combined Table.

Print

You can print the displayed sizing table by clicking Print in the Application Toolbar, then selecting Print.

For More Information about Workflows

Refer to the:

- Peak Scanner™ Software Version 1.0 Reference Guide (P/N 4382253) – Shipped with the software or available online.
- Peak Scanner™ Software Version 1.0 Online Help is accessible from the Peak Scanner™ Software by clicking Help () ► Help Contents in the application toolbar.

© Copyright 2007, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Notice to Purchaser: License Disclaimer

Purchase of this software product alone does not imply any license under any process, instrument or other apparatus, system, composition, reagent or kit rights under patent claims owned or otherwise controlled by Applera Corporation, either expressly or by estoppel.

Applera, Applied Biosystems, and AB (Design) are registered trademarks, and Scanner is a trademark of Applera Corporation or its subsidiaries in the U.S. and certain other countries. All other trademarks are the sole property of their respective owners.

02/2007

www.appliedbiosystems.com



Part Number 4383719 Rev. A